

## EDMUS 4.1 INSTALLATION ON WINDOWS VISTA®

Because of the important differences between Windows Vista and XP, and the number of possible configurations, a specific installation procedure must be followed.

**VERY IMPORTANT:  
READ THESE INSTRUCTIONS CAREFULLY BEFORE YOU INSTALL EDMUS.**

If you are reading this, you are probably in one of these four situations:

- *This is your first installation of EDMUS ever.*  
Please read paragraph (1) below.
- *This is your first installation of EDMUS on Vista, but you have EDMUS on an older computer.*  
Please read paragraph (2) below.
- *You have already installed EDMUS 3 (3.4 to 3.6.1) on Vista, and you want to update it.*  
Please read paragraph (3) below.
- *You have already installed EDMUS 4 (or 3.6.2) on Vista, and you want to update it.*  
Please read paragraph (4) below.

### 1. THIS IS YOUR FIRST INSTALLATION OF EDMUS EVER

***You only have to do a standard installation of EDMUS.***

#### **STEP 1**

Proceed with the installation of EDMUS, using the installer program (this will be a "NEW" installation).

#### **STEP 2**

- Launch EDMUS 4.1. If asked, authorize Omnis Studio to run the studiorg.exe program, this will happen only once.
- Enter your user name and password (by default, "ADMIN" and "ADMIN").
- EDMUS will then automatically import new thesaurus files, this will only take a few seconds.
- Important: enter the EDMUS license information (go to the "Manager" menu -> "License..."), using the EDMUS and Omnis Studio serial numbers that have been emailed to you (the Omnis Studio serial number should begin with "WBRN43").

<b>2. THIS IS YOUR FIRST INSTALLATION OF EDMUS ON VISTA, BUT YOU HAVE EDMUS ON AN OLDER COMPUTER</b>
--

***After installing EDMUS, you will have to transfer your patient data from your older computer to the new one.***

**STEP 1**

Proceed with the installation of EDMUS, using the installer program (this will be a "NEW" installation).

**STEP 2**

Copy your EDMUS datafile (containing your patient data) from the older computer (running Windows XP or 2000) onto the new one.

On your old computer, the datafile is located at:

`C:\Program Files\EDMUS Software\EDMUS 3\Edmus.df1`

or at:

`C:\Program Files\EDMUS Software\EDMUS\Edmus.df1`

On your new computer, the datafile is located at:

`C:\Users\Public\EDMUS Software\EDMUS\Edmus.df1`

Simply replace the Edmus.df1 datafile on your new computer by a copy of the one from your old computer (always keep a backup copy of your older datafile by security).

Note: if the file with a .df1 extension is not named Edmus.df1, rename it Edmus.df1.

**STEP 3**

- Launch EDMUS 4.1. If asked, authorize Omnis Studio to run the studiorg.exe program, this will happen only once.
- EDMUS 4.1 will then convert your old datafile to its new format. This may take a while, depending on the number of patient files you have.
- Enter your user name and password.
- EDMUS will then automatically import new thesaurus files, this will only take a few seconds.
- Important: enter the EDMUS license information (go to the "Manager" menu -> "License.."), using the EDMUS and Omnis Studio serial numbers that have been emailed to you (the Omnis Studio serial number should begin with "WBRN43").

### 3. YOU HAVE ALREADY INSTALLED EDMUS 3 (3.4 to 3.6.1) ON VISTA

***You will have to backup your patient data before installing EDMUS, then transfer these patient data into the new installation.***

#### **STEP 1**

Make a **backup copy** of your old datafile.

Its location will differ according to the configuration of Vista.

- The most likely location is (look there first):

C:\Users\<<your user name>\AppData\Local\VirtualStore\Program Files\  
EDMUS Software\EDMUS 3\Edmus.df1

To go there, you will have to temporarily set Vista to "Show hidden files and folders" (go to "Control Panel" -> "Appearance and Personalization" -> "Folder Options" -> "Show hidden files and folders" -> "View" tab).

- *If not in the first location, your datafile should be at:*

C:\Program Files\EDMUS Software\EDMUS 3\Edmus.df1

But *always* look in the previous location first.

- If Edmus.df1 is present in both locations, the one that contains your patient data is most likely the one in the first location (VirtualStore). Keep a copy of both by security.
- The Edmus.df1 file should be at least 600 kb, usually much more (depending on the number of patient files it contains and on the EDMUS version).

#### **STEP 2**

***(Provided you have found and backed up your EDMUS datafile)***

- Uninstall your old copy of EDMUS 3 (go to "Control Panel" -> "Programs" -> "Uninstall a program" -> "EDMUS 3" -> click on "Uninstall/Change").
- You should also remove the "EDMUS Software" folder that may reside in the VirtualStore folder (put it in the Recycle Bin):

C:\Users\<<your user name>\AppData\Local\VirtualStore\Program Files\  
EDMUS Software

#### **STEP 3**

Launch the installer program again, and this time proceed with the **installation** of EDMUS (select "You have already installed a previous version of EDMUS (3.4 to 3.6.1) on Vista", then "Continue the installation").

#### **STEP 4**

**Replace** the empty Edmus.df1 file located at:

C:\Users\Public\EDMUS Software\EDMUS\Edmus.df1

by a copy of the older one you have just backup up (always keep a copy of the backup file by security).

**STEP 5**

- Launch EDMUS 4.1. If asked, authorize Omnis Studio to run the studiorg.exe program, this will happen only once.
- EDMUS 4.1 will then convert your old datafile to its new format. This may take a while, depending on the number of patient files you have.
- Enter your user name and password.
- EDMUS will then automatically import new thesaurus files, this will only take a few seconds.
- Important: enter the EDMUS license information (go to the "Manager" menu -> "License..."), using the EDMUS and Omnis Studio serial numbers that have been emailed to you (the Omnis Studio serial number should begin with "WBRN43").

**4. YOU HAVE ALREADY INSTALLED EDMUS 4 (OR 3.6.2) ON VISTA**

***You should (as always) make a backup copy of your patient data before updating EDMUS.***

**STEP 1**

Make a **backup copy** of your old datafile.  
Its location is:

C:\Users\Public\EDMUS Software\EDMUS\Edmus.df1

The Edmus.df1 file should be at least 600 kb, usually much more (depending on the number of patient files it contains and on the EDMUS version). Always keep an up-to-date backup copy of this file by security.

**STEP 2**

Launch the installer program again, and this time proceed with the **update** of EDMUS (select "You have already installed EDMUS 4 (or 3.6.2) on Vista", then "Continue the installation", then "UPDATE").

**STEP 3**

- Launch EDMUS 4.1. If asked, authorize Omnis Studio to run the studiorg.exe program, this will happen only once.
- You will then have to enter your Omnis Studio serial number, which has been emailed to you (the Omnis Studio serial number should begin with "WBRN43").
- EDMUS 4.1 will then convert your old datafile to its new format. This may take a while, depending on the number of patient files you have.
- Enter your user name and password.
- EDMUS will then automatically import new thesaurus files, this will only take a few seconds.